



TECHNICAL SPECIFICATIONS

Seat Pad Size: 17.5" wide x 15" deep
(445 mm wide x 381 mm deep)

Width Between Arms: 20" (508 mm)

Footrest Size:
10.5" long x 11" wide
(267 mm long x 279 mm wide)



Stay in the home you love.

Bruno's most popular stairlift, the Elan, can help restore full access to your home. Designed for straight indoor staircases, the Elan offers a mix of high-quality standard functionality, unsurpassed reliability, and ease-of-use design.

With Made-in-the-USA dependability, Bruno helps you stay in your home with comfort and ease.

Learn more about Bruno's best-selling stairlift, the Elan, at www.bruno.com.

ELAN | SRE 3000



Learn more at www.bruno.com



Compact, flip-up design for plenty of open space on stairway.



Two wireless call/sends allow remote control access.



Sensors ensure safety

Sleek appearance





Offset swivel seat makes getting on/off easy.

ELAN FEATURES

- Weight capacity: 300 lb (136 kg)
- Padded, generous-size seat with multiple height adjustments
- Armrest rocker switch for easy operation
- Anodized rail with covered gear rack gives a sleek appearance
- Vertical rail design enables tight installation to within 5" (127 mm) of wall
- Flip up arms, seat, and footrest create plenty of extra space on steps
- Continuous charge strips along rail repower unit wherever it stops
- Two wireless call/sends allow unit to be controlled remotely
- Two 12v batteries ensure dependable performance even in a power outage
- Soft start/stop gives a smooth ride from start to finish
- Obstruction safety sensors, seat belt ensure safety
- Onboard audio diagnostics for easy service
- Gold warranty: five year major components; two year parts

Padded, generous seat



ELAN OPTIONS



Power or manual folding rail for narrow hallway or when doorway is at the bottom of stairs. Manual operation or push button automatic.



Power swivel seat for effortless exit. Armrest switch control or wireless call/send



Power folding footrest automatically flips up/down when seat is raised/lowered. Arm switch control optional.



About Bruno

Bruno Independent Living Aids has helped improve the lives of people with limited mobility for over 30 years. Engineered and manufactured in the USA, Bruno stairlifts and vertical platform lifts help people better access homes, public buildings and businesses, and Bruno's full line of scooter/powerchair lifts help people easily transport their mobility devices. In addition, Bruno's Valet Signature Seating is the world leader in revolutionary turning automotive seats to assist people in getting in and out of their vehicles more easily. A family-owned business, Bruno takes pride in their unwavering focus on quality and safety and their position as industry innovation leader.

Proudly made in the U.S.A.



Bruno Independent Living Aids, Inc.
1780 Executive Drive, Oconomowoc, WI 53066
www.bruno.com

© Bruno Independent Living Aids, Inc. 2015

MKT-SL-1